





Height range 45cm to 98cm

## Medicare 2 Section Couch

Anti Bacterial 240Kg Low Transfer He

A very strong and rigid couch available with a wide range of accessories. Suitable for supporting very heavy patients for examination and minor treatment procedures.



## **Standard Features**

- Safe working load and lifting capacity of 240Kg (528lb
- Gas assisted backrest adjustable from -25° to 75° on
- Choice of hydraulic or electric height elevation
- Retractable wheel design with adjustable foot for stab
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- 63cm wide Easy Clean Upholstery as standard
- Choice of Standard or Standard+ upholstery colours



Low patient transfer height

## **STANDARD** STANDARD+ **Available** Upholstery **Colours:** Black White ——— Light Grey ——— Cream ——— Sky Blue Canard Dark Grev Olive Dark Blue ----- Lotus Green ----- Pepper

	Optional Accessories	
s)	9002	Side support rails (Pair)
	9002P	Side support cushions (Pair)
n	6004	Breathing hole and plug
0	6005	Paper roll holder
	6006	Upper limb support
	6007	Removable head/knee support
ility	6053	Protective base frame cover
	6054	Telescopic IV pole
	6057	Battery back-up facility
	-EPFS	Perimeter Foot Switch for SM2560
	-1WF	Single foot switch for SM2560
	-2WF	Double foot switch for SM2570
	9200/70	70cm Width Easy Clean Upholstery



Negative head section to -25°

