



Clinical Couch Range

Therapy Couch Range



Seating Range



SEERS Medical Ltd Kenton Road, Debenham, Suffolk IP14 6LA, United Kingdom © Copyright SEERS Medical Ltd 2022



info@seersmedical.com www.seersmedical.com

Tel +44 (0)1728 861 590 **Fax** +44 (0)1728 861 868

ISO 13485 ISO 14001

2015 -

- 2016

lssue № 8/22 Part № 06Z-8009-SF

Company Reg. № 6572698 VAT № 931596702 Printed in the United Kingdom

Note¹: The illustrated fabric colours may vary slightly from actual fabric supplied. This occurs due to the four colour printing process. However every attempt has been made to make them as accurate as possible. Our Policy of continuous improvement means that we reserve the right to change and modify designs without prior notice. Note²: Measurements may vary +/- 2° or +/- 5mm due to manufacturing tolerances





Sterling Range

Economy Couches For Patient Examination & Minor Procedures



Height range 46cm to 99cm

Sterling Changing Table

The Sterling Changing Tables provide a variable height support surface to reduce back strain and posture problems for carers and support staff during routine changing procedures in environments such as special schools, specialist adult centres and paediatric units. These tables are available in two sizes and supplied with side support rails as standard. They are fitted with large wheels with brakes for ease of transport and positioning.

Model

SX1047 SX1057 Short table electric - 152cm Long table electric - 184cm

Standard Features

- Safe working load and lifting capacity of 200Kg
- Superb low transfer height for patients
- Choice of 2 lengths, 152cm or 184cm
- 100mm diameter individually locking castors
- Side support rails fitted as standard on all models

- Hand switch fitted as standard on all models
- Standard width upholstery 63cm








Model SX1057 Colour Sky Blue







Height range 46cm to 99cm

Sterling 2 Section Couch

The Sterling 2 Section Couch provides excellent value for money and is suitable for a wide range of examination procedures. The fully adjustable backrest supports the patient for treatment or examinations in the seated position and the backrest can also be adjusted to -20° for neck and head examinations and some manipulative procedures. Additionally, an optional breathing hole and plug can be added for prone treatments. A hand switch is fitted as standard for height control with a foot switch available as an option.



Standard Features

- Safe working load and lifting capacity of 200Kg

- Excellent height range from 46cm to 99cm
- 100mm diameter individual locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- 63cm wide Easy Clean Upholstery as standard
- Choice of Standard upholstery colours



Negative head section to -20°

Available Upholstery Colours:

Colour Dark Blue



Models

SX2007

Model SX2007 Shown with optional breathing hole and plug



Optional Accessories

6004 6003

6005 Paper roll holder Breathing hole and plug Pair of side support rails **9002P** Pair of side support rail cushions -1WF Single foot switch operation **6007** Matching headrest/ knee support




Low patient transfer height

Dark Grev

5





Sterling 3 Section Couch

Like the 2 Section Model, the Sterling 3 Section Couch is excellent value for money and the 3 section format provides a more versatile and comfortable position for seated procedures and examinations. The backrest support can be raised and fully adjusted from horizontal to +85° and the foot section is adjustable from -20° to +40° and can also be used as the head end for prone treatments with the option of adding a breathing hole and plug. The couch is therefore suitable for a wide range of examination and minor treatment procedures. A hand switch is fitted as standard for height control with a foot switch available as an option.




- Safe working load and lifting capacity of 200Kg

- 100mm diameter individual locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- 63cm Wide Easy Clean Upholstery as standard
- Choice of Standard upholstery colours







Optional Accessories

6005 Paper roll holder 6004 Breathing hole and plug Pair of side support rails 6003 **9002P** Pair of side support rail cushions -1WF Single foot switch operation **6007** Matching headrest/ knee support

7





Optional Accessories

6003

-1WF

Fold-down side support rails for 2 & 3 section couches, providing patient support

Single foot switch option

(alternative to standard single hand switch).

6007



easily attached to the couch support for lower limbs.



6004

Breathing hole & plug for improved patient comfort for

9002P

rails for additional comfort, support and protection.



6005

The paper roll holder upholstered surface.

Mundial Standard Range

A range of globally popular, hard-wearing vinyl colours that will retain performance for many years. Abrasion Resistance to 150,000 cycles. Available as standard on all upholstered couches and seating products.





Sky Blue

Canard



Dark Grey

If you have any queries on any of the optional accessories listed above please contact us for further details. Please note some of the accessories are supplied as standard with specific couches so please ensure you check before purchasing your chosen product as all other options come at an additional cost.

Please Note: The illustrated vinyl colours may vary slightly from actual vinyl supplied. This occurs due to the four colour printing process. However every attempt has been made to make them as accurate as possible.

Dark Blue



Black