











## Therapy Deluxe Couch



A very strong and rigid couch with full length top frame for added stability. These models are available with a choice of three head section designs with the option of split leg facility and have the perimeter foot switch operation supplied as standard.

## Models

Standard head design ST3347

ST3347S

ST3348 Basic head design

ST3348S

ST3349 Plus head design

ST3349S Plus head design with split leg



## **Standard Features**

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric height operation via integral perimeter foot switch
- Excellent height range from 46cm to 99cm

- Gas assisted backrest and head sections incorporate anti-trap design
- Choice of standard, basic or plus head design with adjustable armrests
- Choice of standard or split leg foot section
- Breathing hole and plug fitted as standard

## **Optional Accessories**

**6007** Removable Head/Knee Support **6076** Memory Foam Upholstery



**Available Upholstery Colours:** 

