

Easy Clean Upholstery (ECU)

Fitted as Standard on





Sterling 2 Section Couch



Sterling 3 Section Couch



Therapy 2 Section Couch - Standard Head ONLY



Medicare 2 Section Couch



Therapy 3 Section Couch - Standard Head ONLY



Medicare 3 Section Couch

70cm Wide Easy Clean Upholstery



All CLINNOVA Couches

70cm Wide Easy Clean Upholstery Available at an Additional Cost on:

- Medicare 2 Section Couches
- Medicare 3 Section Couches
- Therapy 2 Section Couches
 - Standard & Basic Head ONLY
- Therapy 3 Section Couches









Benefits

- No Visible Staples reduces risk of injury whilst cleaning
- No exposed staples or folded material, means bacteria can no longer grow
- The vacuum formed upholstery finish is made from strong recycled ABS material. Making it highly durable and long lasting
- Easier surface to clean and maintain
- 90° return prevents liquid
- Designs registered with the Intellectual Property Office 6038182, 6038183, 6038184 & 6038185
- Secured with super strong, flush-fitting poppers

