


## Easy Clean Side Cushions

SEERS Medical have innovatively designed these highly durable side support cushions with a staple-free easy clean finish. The unique backing prevents bacteria from growing in hard-to-reach spots, creating an easier surface to clean.

### Models

**NV9002P** Easy Clean Side Cushions

*\*Please check prior to ordering NV9002P that the couch in question is being fitted with or is already fitted with one of the Side Support Rail Models listed.*



Back view of NV9002P

No Visible Staples

### Standard Features

- Provides a smooth, quick-wipe surface for easy cleaning
- No visible staples - reduces risk of injury whilst cleaning
- 90° return prevents liquid ingress
- Secured with super strong, flush-fitting poppers
- Highly durable and rigid
- Ultraviolet resistant to sunlight to retain colour
- Vacuum formed for reliable quality & finish
- Manufactured in the UK from recycled ABS, a strong material commonly used in the automotive industry
- Perfectly matched with anti-bacterial vinyl for the ultimate infection control solution

### Order Now

NV9002P is compatible with **new & existing** SEERS Couches that are being supplied with Side Support Rails:

### Side Support Rail Models

<b>NV9002</b>	Side Support Rails for Easy Clean
<b>NV9003</b>	Side Support Rails for Easy Clean
<b>ST6002</b>	Side Support Rails
<b>ST6003</b>	Side Support Rails
<b>ST6073</b>	Side Support Rails for Bariatric 3 Section
<b>ST6074</b>	Side Support Rails for Mat/Bobath/Bariatric 2 Section



Front view of NV9002P  
Colour: Sky Blue



Side Support Rails,  
Required for use with NV9002P

## ST6002P / ST6003P:

### Backing Material:

(Hard-wearing, waterproof vinyl, with an abrasion resistance of 150,000 cycles)  
Anti-Bacterial Vinyl

### Vinyl Backing is Secured with Staples:

Bacteria could potentially grow under the staples over time & if a staple becomes loose this could create a risk when cleaning



Crib 5

### ST6002P:

Compatible with only ST6002

### ST6003P:

Compatible with only ST6003




## NV9002P:

### Easy Clean Backing is Secured with Super Strong Flush-Fitting Poppers:

Providing a smooth quick-wipe surface for easy cleaning

### Backing Material:

(Made from Recycled ABS Material, highly durable and will last for years)  
Easy Clean Backing



### NV9002P:

Compatible with ST6002, ST6003, NV9002, NV9003, ST6073 & ST6074



Crib 5