





CLINNOVA Gynae 2









The CLINNOVA Gynae 2 is supplied with Side Support Loops, fully adjustable Leg Supports and a crescent shaped seat design for colposcopy & gynaecology procedures. Supplied with a choice of wheel systems, with a classic or premium base.



Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)

- Variable height range of 46cm to 95cm for individual working heights
- Choice of independent, retractable or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Includes a hand switch as standard
- A pair of multi-adjustable leg supports and arm loops supplied as
- Seat section lifts for easy cleaning of framework
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Paper Roll Holder Fitted as Standard

Optional Accessories

6207 Double foot switch

Removable head/knee support

Removable leg extension

Battery back-up facility

Telescopic IV pole 6054



Colour Black

36



CLINNOVA Gynae 2 - Models Available



NV8473-CLS

Electric Height Elevation Electric Backrest Independently Braking Castors Classic Base



NV8573-CLS

Electric Height Elevation Electric Backrest Retractable Castors Classic Base



NV8673-CLS

Electric Height Elevation Electric Backrest Central Locking Castors Classic Base



NV8473-PRM

Electric Height Elevation Electric Backrest Independently Braking Castors Premium Base



NV8573-PRM

Electric Height Elevation Electric Backrest Retractable Castors Premium Base



NV8673-PRM

Electric Height Elevation Electric Backrest Central Locking Castors Premium Base

38