



CLINNOVA Therapy Hygiene



A selection of sturdy, variable height changing tables suitable for individuals with additional living needs. Popular within schools for children with special educational needs, but also suitable for adolescents and adults.



Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Choice of Small, Medium or Large upholstery surfaces
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Lock & Go transport mode for easy moving
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy

- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples



Side Support Rails Fitted as Standard



46cm to 95cm

Optional Accessories

6200 Single foot switch (*NV1561 only*) **9002P** Side support rail cushions (Pair)



Model NV1551C-CLS Colour Sky Blue



CLINNOVA Therapy Hygiene - Models Available



NV1551A-CLS

Hydraulic Height Elevation Large Classic Base



NV1561A-CLS

Electric Height Elevation Large Classic Base



NV1551B-CLS

Hydraulic Height Elevation Medium Classic Base



NV1561B-CLS

Electric Height Elevation Medium Classic Base



NV1551C-CLS

Hydraulic Height Elevation Small Classic Base



NV1551A-PRM

Hydraulic Height Elevation Large Premium Base



NV1561A-PRM

Electric Height Elevation Large Premium Base



NV1551B-PRM

Hydraulic Height Elevation Medium Premium Base



NV1561B-PRM

Electric Height Elevation Medium Premium Base



NV1551C-PRM

Hydraulic Height Elevation Small Premium Base





NV1561C-CLS

Electric Height Elevation Small Classic Base



NV1561C-PRM

Electric Height Elevation Small Premium Base