









## **CLINNOVA** Versa 2

Models





**Bacterial** 



The CLINNOVA Versa 2 is designed for podiatry, chiropody & lower limb procedures. The individually adjustable footrests can be extended to accommodate for taller patients. A pair of fold-down armrests are fitted for patient security and comfort.



# **Standard Features**

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)

- Anti-trap gas assisted sections prevent risk of crushing and harm to patient, user or furniture
- Electric height elevation, raises from 46cm to 95cm in 15 seconds
- Choice of Retractable, independent or central locking castors
- A pair of side support loops supplied as standard for patient security
- and increase product life expectancy
- Frame design enables use of mobile hoist
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

# **Optional Accessories**

6004	Breathing hole and plug
6207	Double foot switch
6006	Upper limb support
6041	Side support loops
6044	Lithotomy Stirrups
6045	Gynaecology Leg supports
6007	Removable head/knee supr

**6057** Battery back-up facility Telescopic IV pole



Model NV0675-CLS Colour Light Grey



Paper Roll Holder

Fitted as Standard



# **CLINNOVA** Versa 2 - Models Available



## **NV0475-CLS**

Electric Height Elevation Electric Head Section Gas Assisted Foot Section Independently Braking Castors Classic Base



# **NV0575-CLS**

Electric Height Elevation Electric Head Section Gas Assisted Foot Section Retractable Castors Classic Base



## **NV0675-CLS**

Electric Height Elevation Electric Head Section Gas Assisted Foot Section Central Locking Castors Classic Base



#### **NV0475-PRM**

Electric Height Elevation Electric Head Section Gas Assisted Foot Section Independently Braking Castors Premium Base



# **NV0575-PRM**

Electric Height Elevation Electric Head Section Gas Assisted Foot Section Retractable Castors Premium Base



## **NV0675-PRM**

Electric Height Elevation Electric Head Section Gas Assisted Foot Section Central Locking Castors Premium Base

50 51