





CLINNOVA Mobile 3 Section







A very strong and rigid couch, fitted with large castors which make movement much easier. Fitted with a paper roll holder and side support rails as standard. Suitable for heavy patients for examination and minor treatment procedures.



Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)

- Anti-trap gas assisted sections prevent risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Choice of 125mm independent or central locking castors
- and increase product life expectancy
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

6004 Breathing hole and plug

Single foot switch (NV3X60 only)

Double foot switch (NV3X70 only)

9003P Side support rail cushions (Pair)

6006 Upper limb support

Removable head/knee support

Battery back-up facility

6054 Telescopic IV pole



Paper Roll Holder Fitted as Standard



Colour Sky Blue

16 17



CLINNOVA Mobile 3 Section - Models Available



NV3450-CLS

Hydraulic Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Independently Braking Castors Classic Base



NV3450-PRM

Hydraulic Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Independently Braking Castors Premium Base



NV3460-CLS

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Independently Braking Castors Classic Base



NV3460-PRM

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Independently Braking Castors Premium Base



NV3470-CLS

Electric Height Elevation Electric Head Section Electric Foot Section Independently Braking Castors Classic Base



NV3470-PRM

Electric Height Elevation Electric Head Section Electric Foot Section Independently Braking Castors Premium Base



NV3650-CLS

Hydraulic Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Central Locking Castors Classic Base



NV3650-PRM

Hydraulic Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Central Locking Castors Premium Base



NV3660-CLS

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Central Locking Castors Classic Base



NV3660-PRM

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Central Locking Castors Premium Base



NV3670-CLS

Electric Height Elevation Electric Head Section Electric Foot Section Central Locking Castors Classic Base



NV3670-PRM

Electric Height Elevation Electric Head Section Electric Foot Section Central Locking Castors Premium Base

18 19