



CLINNOVA Therapy 2 Section (Basic)



A very strong and rigid couch available with a wide range of accessories. Suitable for supporting very heavy patients for examination and minor treatment procedures.

Models

See Page 10-11 for available models



Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted head section, adjustable from -25° to +45°
- Anti-trap gas assisted head section prevents risk of crushing harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working he
- Linked retractable wheels with adjustable foot for uneven flo
- Lock & Go transport mode for easy moving
- Purpose-made precision bearings ensure quiet, reliable ope and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours



Offset gas strut for unobstructed patient view through breathing hole



	Opt 9001 6200 9002	ional Accessories Perimeter foot switch (NV2568 only) Single foot switch (NV2568 only) Side support rails (Pair)
and	9002P 9005 6007	Side support rail cushions (Pair) Paper Roll Holder Removable head/knee support
eights oors	9008	Integral Face Cradle Cushion
ration		
of		



Model NV2558-CLS Colour Burgundy



CLINNOVA Therapy 2 Section (Basic) - Models Available



NV2558-CLS

Hydraulic Height Elevation Gas Assisted Head Section Classic Base



NV2568-CLS

Electric Height Elevation Gas Assisted Backrest Hand or Foot Switch Classic Base





NV2558-PRM

Hydraulic Height Elevation Gas Assisted Head Section Premium Base



NV2568-PRM

Electric Height Elevation Gas Assisted Head Section Hand or Foot Switch Premium Base



NV2568-CLS-EPFS

Electric Height Elevation Gas Assisted Head Section Perimeter Foot Switch Classic Base



NV2568-PRM-EPFS

Electric Height Elevation Electric Head Section Perimeter Foot Switch Premium Base