





CLINNOVA Therapy 3 Section (Plus)







A very strong and rigid couch available with a wide range of accessories. The additional section provides extra versatility and allows a seated position. Suitable for supporting very heavy patients.



Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted head section, adjustable from -25° to +45°
- Gas assisted foot section, adjustable from 0° to +85°
- Anti-trap gas assisted head section prevents risk of crushing
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working height
- Linked retractable wheels with adjustable foot for uneven floors
- Lock & Go transport mode for easy moving
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hois
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colour

Optional Accessories


9001 Perimeter foot switch (NV3569 only)6200 Single foot switch (NV3569 only)

9003 Side support rails (Pair)

9003P Side support rail cushions (Pair)

9005 Paper Roll Holder

6007 Removable head/knee support9008 Integral Face Cradle Cushion



Adjustable armrests included as standard on Plus models



Model NV3559-CLS Colour Burgundy

24 25



CLINNOVA Therapy 3 Section (Plus) - Models Available



NV3559-CLS

Hydraulic Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Classic Base



NV3569-CLS

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Hand or Foot Switch Classic Base



NV3569-CLS-EPFS

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Perimeter Foot Switch Classic Base



NV3559-PRM

Hydraulic Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Premium Base



NV3569-PRM

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Hand or Foot Switch Premium Base



NV3569-PRM-EPFS

Electric Height Elevation Gas Assisted Head Section Gas Assisted Foot Section Perimeter Foot Switch Premium Base

26 27