



Medicare Bariatric 2 Section Couch



The Bariatric 2 Section Couches feature electric height and backrest operation for effortless patient positioning for various medical examinations and treatment procedures. There is a choice of widths available - Models SM2573 and SM2673 are available with a width of 80cm and Models SM4542 and SM4642 which have a wide frame, are available at 105cm as standard. All couches provide a safe working load and lifting capacity of 325Kg (715lbs).

Standard Features

- Safe working load and lifting capacity of 325Kg (715lbs)
- Choice of retractable wheels or central locking castors
- Electric backrest adjustable from -15° to +89°
- Electric height elevation from 45cm to 98cm
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Choice of 80cm or 105cm width upholstery
- Hand switch fitted as standard
- Choice of Standard and Standard+ upholstery colours

Optional Accessories

- 6004** Breathing hole and plug
- 6005** Paper roll holder
- 6074** Side support rail (Pair)
- 9002P** Side support rail cushions (Pair)
- 6076** Memory foam upholstery (60mm)
- 2WF** Double foot switch

Models

- SM2573** Electric with electric backrest 80cm wide. Retractable wheel system
- SM2673** Electric with electric backrest 80cm wide. Central Locking Castors
- SM4542** Electric with electric backrest 105cm wide. Retractable wheel system
- SM4642** Electric with electric backrest 105cm wide. Central Locking Castors



Model SM4542
Shown with 6074 Side Rails and 9002P Side Cushions



Model SM2573
Colour Dark Blue



Model SM4642
Colour Sky Blue

Available Upholstery Colours:

