





Standard Dimensions

\*When fitted with 9300/70 Easy Clean Upholstery

## Medicare 3 Section Couch









The Medicare 3 Section Couch provides a very stable and versatile support surface for prone, supine and seated treatment procedures and examinations. Both the elevating sections can be set to a wide range of angles for specific patient positioning.

## Models

SM3550 Hydraulic with gas assisted backrest SM3560 Electric with gas assisted backrest SM3570



## **Standard Features**

- Safe working load and lifting capacity of 240Kg (528lbs)

- Choice of hydraulic or electric height elevation

- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- 63cm wide Easy Clean Upholstery as standard
- Choice of Standard or Standard+ upholstery colours

## **Optional Accessories**

**9003** Side support rails (Pair)

**9002P** Side support cushions (Pair)

Breathing hole and plug

Paper roll holder

Upper limb support

Removable head/knee support

Protective base frame cover

Telescopic IV pole

**6057** Battery back-up facility

Perimeter Foot Switch for SM3560

Single foot switch for SM3560

Double foot switch for SM3570

9300/70 70cm Width Easy Clean Upholstery







**Available Upholstery Colours:** 

