





Medicare Echocardiography Couch









Designed specifically for Echocardiography, this couch has electric operation of height, tilt and backrest angle. The backrest has a cut-away design for improved application of the transducer during scanning procedures and a specially designed stool for the clinician is also provided. The couch is also supplied with an integral footrest for the clinician, as well as a patient hand-grip and head cushion for the patients comfort and security.

## Models

SM3590

Electric height, tilt and backrest. Electric height, tilt and backrest. Central locking wheel system

SM3690



## **Standard Features**

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric height, tilt and backrest adjustment

- Electric trendelenburg tilt to -15°
- incorporating anti-trap design
- Ergonomic backrest designed for improved transducer positioning
- Integral footrest for clinician when traditional right-hand scanning, together with specially designed stool
- Hand-grip and side support with cushion for patient positioning
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Paper Roll Holder and Head Support Cushion fitted as standard
- Retractable wheel design with adjustable foot for stability on
- Choice of Standard and Standard+ upholstery colours

## **Optional Accessories**

**-3WF** Triple foot switch operation Battery back-up facility

**6123** Fully adjustable arm support for

the clinician



Special design stool provided for clinician, included as standard





Model 6123 optional adjustable arm support for clinician

**Available Upholstery Colours:** 

