





Deluxe Saddle Chairs









The Deluxe Saddle Chairs are commonly used in Sonography, Dentistry and Surgical environments. The ergonomic design promotes positive posture by encouraging pelvic tilt whilst the user straddles the stool with both feet firmly on the floor and body weight is distributed through the thighs. The 2 lever mechanism operates and adjusts the seat height and tilt. All models feature a polished aluminium five star base and castors, with option of foot height control.

Models

MC6107

Deluxe saddle chair - standard model



High quality upholstery finish



Model MC6107 Colour Dark Blue

Available Upholstery Colours:

