







## **Round Medical Stools**







The clinical stool range has been designed to offer the clinician a variety of seat options, promoting optimum posture during various procedures. The stools are suitable for use in a wide range of clinical environments. The standard model provides a height range of 45cm to 59cm and the high model provides a height range of 54cm to 74cm. High models are supplied with an adjustable foot support ring as standard. The foot height control options provides hands free height adjustment.

## Models

MC6164 MC6165 Round medical stool - standard model Round medical stool - high model Foot support ring fitted as standard.



All round height adjustment ring, fully upholstered under side for ease of cleaning



Model MC6164
Colour Diabolo Mandarine
Colour available at additional cost



Model MC6165 with 6174 Foot Height Control
Colour Diabolo Mandarine

## **Available Upholstery Colours:**

