





Sterling Couch Range



Therapy Couch Range



Clinical Couch Range





LRQ00001104

SEERS MEDICAL LtdKenton Road,
Debenham, Suffolk

IP14 6LA, United Kingdom © Copyright SEERS Medical Ltd 2022

Company Reg. № 6572698 VAT № 931596702 Printed in the United Kingdom



Tel +44 (0)1728 861 590 **Fax** +44 (0)1728 861 868

info@seersmedical.com www.seersmedical.com

> Issue № 8/22 Part № 06Z-8008-SF

Note¹: The illustrated fabric colours may vary slightly from actual fabric supplied. This occurs due to the four colour printing process. However every attempt has been made to make them as accurate as possible. Our Policy of continuous improvement means that we reserve the right to change and modify designs without prior notice. Note²: Measurements may vary $+/-2^{\circ}$ or +/-5mm due to manufacturing tolerances



Seating Range

Clinical & Rehabilitation Seating Products





Introduction

Since **SEERS** Medical was founded in 2008, the company has always offered some basic chairs and stools to complement its range of medical furniture. Due to the increasing demand from customers to provide a complete range of clinical and patient seating, **SEERS** has developed a dedicated medical seating range to meet these growing requirements. In this brochure **SEERS** presents many new designs as well as a number of other complementary medical seating products that we distribute, providing a varied choice of models and options for various medical procedures and examinations.

The development of these products has been achieved by applying our many years of experience and knowledge in this market. Much thought has been given to the design, to benefit both the comfort and security of the patients as well as providing manual handling solutions which promote good posture and ergonomics to all clinical staff when they are using this equipment. With the management and control of our own assembly line at our premises in Debenham, Suffolk in the UK, our couches comply to the CE medical devices directive 93/42/EEC, ensuring we supply high quality products every time, with short lead times to meet the tightest of deadlines. We trust you will experience this yourself when you purchase a **SEERS** Medical product.

SEERS Medical are ISO 14001 and ISO 13485:2016 accredited.

Contents

2-3	Introduction & Technical Information
4-5	Proximity & Keyhole Shaped Stools
6-7	Dual Curve Stools
8-9	Round & Square Medical Stools
10-11	Deluxe Saddle Chairs & Saddle Stools
12-13	Economy Saddle Stool & Sonographers Chair
14-15	Surgeons & Sonographer Chairs
16-17	Operators Stool & Contoured Medical Chair
18-19	Operators Chairs
20-21	General Medical Chairs
22-23	Laboratory Chairs & Stools
24-25	Multi Procedures Chair
26-27	Ophthalmology & Phlebotomy Chair
28-29	Optional Accessories & Vinyl Colour Options



Anti-microbial



Crib 5 Fire retardant upholstery



Complies with EC Directive 93/42/EEC



Intelligent Ergonomic Designs

Technical Information

The **SEERS** Medical Seating range includes a wide variety of ergonomically designed chairs and stools designed for clinicians, surgeons, therapists and more. All **SEERS** Seats are manufactured in Debenham, Suffolk and are individually inspected before being packed. This range complements the existing **SEERS** Clinical & Therapy Ranges, providing comfort and support to the person performing procedures. The collection of seats caters for a large variety of healthcare environments, including hospitals, general practices and therapy clinics.







Here's what some of our customers think...

$\triangle \triangle \triangle \triangle \triangle \triangle$

"Excellent service & great products"

– Sandra

$\triangle \triangle \triangle \triangle \triangle \triangle$

"High quality & durable equipment"

– Cynthia

$\triangle \triangle \triangle \triangle \triangle \triangle$

"No backache since purchasing one of your chairs"


- Michael

"Great products with a reliable and friendly service"

- Martina











Proximity Stools







Keyhole Shaped Stool







This stool is intelligently designed to have one flat side, allowing the clinician to get in close proximity to the patient during a procedure. The stool is available in both standard and chrome bases, both with an excellent height range of 55-81cm. Included as standard with **SEERS** Echocardiography Couches.

Models

MC6105

Proximity Stool with black standard base and glides





An intelligently & ergonomically designed stool for use in a multitude of healthcare applications - from dental practices to tattoo studios. With a 150kg safe working load, the stool is suitable for a range of users and is guaranteed for 5 years.

Models

ST6109

Keyhole Shaped Stool - Standard Model





Model ST6109 Colour Cream

Available Upholstery Colours:























Economy Dual Curve Stool







Ergonomic Dual Curve Stools









The dual curve medical stools feature an ergonomically designed seat pad, the curvature and static tilt of the seat encourages pelvic tilt, prevents slouching and promotes positive posture. An ideal suturing stool for clinical environments such as Delivery Suites, Theatres and A&E Departments, where the clinician is required to work in close proximity to the patient during these procedures.

Models

ST6168

Dual Curve Medical Stool





The dual curve medical stools feature an ergonomically designed seat pad, the curvature and static tilt of the seat encourages pelvic tilt, prevents slouching and promotes positive posture. An ideal suturing stool for clinical environments such as Delivery Suites, Theatres and A&E Departments, where the clinician is required to work in close proximity to the patient during these procedures. The foot height control option provides hands free height adjustment during these procedures.

Models

MC6168 MC6169 Dual curve medical stool - standard model Dual curve medical stool - high model Foot support ring fitted as standard.









Colour Diabolo Tilleul



Model MC6169 with 6174 Foot Height Control Colour Diabolo Tilleul


Available Upholstery Colours:





















Round Medical Stools







Square Medical Stools







The clinical stool range has been designed to offer the clinician a variety of seat options, promoting optimum posture during various procedures. The stools are suitable for use in a wide range of clinical environments. The standard model provides a height range of 45cm to 59cm and the high model provides a height range of 54cm to 74cm. High models are supplied with an adjustable foot support ring as standard. The foot height control options provides hands free height adjustment.

Models

MC6164 MC6165

Round medical stool - standard model Round medical stool - high model Foot support ring fitted as standard.



All round height adjustment ring, fully upholstered under side for ease of cleaning



Model MC6164 Colour Diabolo Mandarine



Model MC6165 with 6174 Foot Height Control Colour Diabolo Mandarine

These models provide the clinician with the option of a square seat pad promoting positive posture by allowing the legs to fall on either side of the seat corners, thus enabling the clinician to work in close proximity to the patient during an examination or procedure.

Models

MC6166 MC6167 Square medical stool - standard model Square medical stool - high model Foot support ring fitted as standard.






Model MC6166 Colour Diabolo Oriflamme



Foot height control option available


Available Upholstery Colours:




















Deluxe Saddle Chairs









Ergonomic Saddle Stools



The ergonomic saddle stools have a similar design to the deluxe saddle chairs, designed to encourage positive posture

throughout seated procedures. The 2 lever mechanism operation provides independent adjustment for the seat height







The Deluxe Saddle Chairs are commonly used in Sonography, Dentistry and Surgical environments. The ergonomic design promotes positive posture by encouraging pelvic tilt whilst the user straddles the stool with both feet firmly on the floor and body weight is distributed through the thighs. The 2 lever mechanism operates and adjusts the seat height and tilt. All models feature a polished aluminium five star base and castors, with option of foot height control.

Models

MC6107

Deluxe saddle chair - standard model



High quality upholstery finish



Model MC6107 Colour Dark Blue

Models

MC6144 MC6145 Foot support ring fitted as standard.


and tilt. All models feature a polished aluminium five star base and castors.





Colour Lotus Green

Available Upholstery Colours:









ST6144 45cm to 59cm











Economy Saddle Stool







Economy Sonographers Chair







The economy saddle stool ensures comfortable posture for therapists and clinicians during procedures. The curved seat pad design provides a similar ergonomic seated position to the general saddle stool. The lever mechanism operation provides adjustment for the seat height. All models feature a black plastic base, gas lift and castors.

Models

ST6144

Economy Saddle Stool





The economy sonographers chair features an ergonomically designed seat pad. The chair has a safe working load of 150kg and comes with a 5 year warranty. The height and backrest are adjustable via the levers under the seat.

Models

ST6140


Economy Sonographers Chair







Available Upholstery Colours:

















Surgeons & Sonographers Chair









The surgeons and sonographers chairs feature an ergonomically designed seat pad and smaller backrest that will not impede the movement of the arm during scans and procedures. The curvature design and seat pad adjustment encourages pelvic tilt, prevents slouching and promotes positive posture. The chair is standard with a small, compact base to reduce trip hazards and features a wipe clean upholstered underside.

Models

MC6158 MC6159 Surgeons chair - standard model Foot support ring fitted as standard



Fully upholstered backrest and seat pad for ease of cleaning



Colour Diabolo Oriflamme

Features

- Wide choice of models available
- Ergonomically designed seat
- Compact backrest for improved scanning procedures
- Weight capacity and safe working load of 150Kg
- Anti-bacterial and waterproof vinyl upholstery
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- **6178** Anti-static castors
- **6455** Load locking castors
- **6172** Manual locking castors
- Glides for added stability on high models
- Height adjustable armrests
- **6174** Foot height control operation
- **6176** Foot support ring for standard model
























Operators Stools







A general operators stool with anti-microbial vinyl, ideal for a multitude of medical and clinical environments. Available in both Standard & High models. The high model features a foot support ring as standard.

Models

ST6103 ST6104

Standard Operators Stool
High Operators Stool with Foot Support

Ring & Glides





New Easy Clean Upholstery System to maximise infection prevention



Contoured Medical Chair









A versatile chair with a contoured seat pad and raised front to promote positive posture whilst seated. The 3 lever mechanism provides independent positioning of the backrest angle, seat pad tilt and chair height. The contoured backrest is ergonomically designed to provide lumbar support for the operator and has the option of foot height control.

Models

MC7212 MC7213 Contoured medical chair - standard mode Contoured medical chair - high model Foot support ring fitted as standard







Available Upholstery Colours:



Available Upholstery Colours:




Canard Dark Blue Dark Grey Black U White Light Grey Cream Olive Lotus Green Pepper Brown Poterie Red Burgundy Sky Blue Canard Dark Grey Black U White Light Grey Cream Olive Lotus Green Pepper Brown Poterie Red Burgundy T















49

Operators Chairs







A general operators chair with anti-microbial vinyl, ideal for a multitude of medical and clinical environments. Included with Medicare Deluxe Gynaecology Couch & various Clinical Packages.

Models

ST6100 ST6101 Standard Operators Chair High Operators Chair with Foot Support Ring & Glides





Premium Operators Chair







The premium operators chair features anti-microbial vinyl, ideal for a wide range of clinical environments. The 3 Lever Mechanism allows the user to find their desired height, backrest angle and pelvic tilt, encouraging posture and ensuring comfort.

Models

ST6150 Premium Operators Chair - StandardST6151 Premium Operators Chair - High







Available Upholstery Colours:




Available Upholstery Colours:














General Medical Chairs









Suitable for a wide range of medical procedures, these chairs provide comfort and support for clinicians sitting for long periods. The deep padded round seat features a pressure reducing foam membrane and in combination with the fully adjustable backrest, promotes good posture and comfort. There is a choice of standard or high models, with high models also featuring a Foot Support Ring as standard.

Models

MC6178 MC6179

20

Round medical chair - standard model Round medical chair - high model Foot support ring fitted as standard





Features

- Deep padded seat with pressure reducing foam
- Suitable for long periods of sitting
- Weight capacity and safe working load of 150Kg
- Wide choice of models available
- Anti-bacterial and waterproof vinyl upholster
- Restraint castors fitted as standard
- Choice of Standard or Standard+ upholstery colour

Optional Accessories

- **6178** Anti-static castors
- **6455** Load locking castors
- **6172** Manual locking castors
- Glides for added stability on high models
- **6175** Height adjustable armrests
- 6174 Foot height control operation
- **6176** Foot support ring for standard model







MC6112 43cm to 57cm

MC6115 52cm to 72cm



MC6111 43cm to 57cm

MC6114 52cm to 72cm



MC6113 41cm to 55cm

MC6116 52cm to 72cm



Laboratory Chairs & Stools







Designed for use within laboratories or clean room environments. Backrests and seat pads are all manufactured using hard wearing chemical resistant polyurethane. All chairs and stools are height adjustable and can be fitted with optional foot height control and a wide choice of castors.

Models

MC6111	Lab chair - standard model
MC6112	Lab stool - standard model
MC6113	Lab saddle stool - standard model
MC6114	Lab chair with foot support ring
MC6115	Lab stool with foot support ring
MC6116	Lab saddle stool with foot support ring





Model MC6114

Features

- Designed for laboratories and cold rooms
- Hard wearing chemical resistant polyurethane (PU)

Optional Accessories

- **6178** Anti-static castors
- **6455** Load locking castors
- **6172** Manual locking castors
- Glides for added stability on high models
- Height adjustable armrests
- Foot height control operation
- **6176** Foot support ring for standard model



Model MC6113




Model MC6115









Multi Procedures Stool & Chair









Ideally suited for Operating Suites, Ultrasound, Dentistry and Ophthalmic Departments. The 360° swing around arm provides the clinician with an armrest or torso support, using a single movement and offering additional stability during procedures. The ratchet system allows the armrest to be adjusted towards or away from the clinician, providing the user with a variety of options for different procedures. The height of the arm can also be adjusted, ideal if it is being used to support the torso. The backrest models provide additional support and comfort. All models are available with a foot height control as an option.

Models

MC6160 Chair with backrest - standard model
MC6161 Chair with backrest - high model
Foot support ring fitted as standard

MC6162 Stool - standard mode MC6163 Stool - high model

Foot support ring fitted as standard



Fully adjustable 360° swing around arm fitted as standard on all models



Features

- Wide Choice of Models Available
- Intelligent & Ergonomic Design
- 360° Swing around arm fitted on all models and can be used as a backrest or armrest
- Safe Working Load of 150kg
- Upholstery vinyl is Anti-Bacterial, Fireproof to Crib 5 and Waterproof
- Glides (Static Feet) fitted as standard
- Available in a wide range of colours
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- **6178** Anti-static castors
- **6455** Load locking castors
- **6172** Manual locking castors
- 173 Glides for added stability on high models
- **6174** Foot height control operation
- **6176** Foot support ring for standard model






Foot Support Ring fitted as standard on Models MC6161 & MC6163

Model MC6161 Colour Diabolo Cobalt Colour available at additional cost


















Ophthalmology Chair









The chair has a low seat height for patient mount and dismount and the hydraulic height adjustment enables the clinician to alter the position of the seated patient to the desired height before carrying out procedures and examinations. The arms provide maximum support for the patient as they mount and dismount the chair. The headrest is able to adjust in height and angle for optimum comfort and positioning.

Models

MC6172

Ophthalmology Chair - Hydraulic



Fully adjustable headrest



Phlebotomy Chair









The hydraulic height adjustable phlebotomy chair enables the operator to easily and accurately position the patient for blood sampling procedures. The chair is a compact and space saving alternative to electrically operated chairs and couches. Supplied with a sturdy base on glides (static feet) to ensure the patient is secure. Integral arm supports are fitted as standard to assist the patient in mounting and dismounting the chair. The chair has a height range of 50cm to 64cm. These models are available with a choice of single or dual phlebotomy armrests.

Models

MC6154

MC6155

Hydraulic phlebotomy chair with single adjustable armrest. Hydraulic phlebotomy chair



Fully adjustable moulded arm supports



Available Upholstery Colours:















Optional Accessories



6178

A set of anti-static castors.



6455

Load locking castors automatically lock when weight is applied to the chair or stool.



6172

Manual locking castors provide the user with the option to manually lock each castor.



6173

Glides (also know as static feet) enable the chair or stool to be kept in a fixed position at all times



6175

with black swivel pads provide comfort for the user.



6348

flexibility and enables the user to fold the arms down to 90° for convenience.



6174

Foot height control (FHC) enables the height to be altered via foot operations as opposed to hand operation



6176

Foot support ring (Fitted as standard on high models). Not to be used with Model 6174 Foot height control.

If you have any queries on any of the optional accessories listed above please contact us for further details. Please note some of the accessories are supplied as standard with specific chairs so please ensure you check before purchasing your chosen product as all other options come at an additional cost.

Mundial Standard Range

A range of globally popular, hard-wearing vinyl colours that will retain performance for many years. Abrasion Resistance to 150,000 cycles. Available as standard on all upholstered couches and seating products.



Mundial Standard+ Range*

A selection of unique coloured, hard-wearing vinyl colours that will retain performance for many years. Abrasion Resistance to 150,000 cycles.

* Available at an additional cost for all Ergonomic Seating models, excluding Lab/PU Seating.



Other colours are available on request

Please Note: The illustrated vinyl colours may vary slightly from actual vinyl supplied. This occurs due to the four colour printing process, however every attempt has been made to make them as accurate as possible.