





Therapy 3 Section Plus Head Couch



Designed with a short head section for head and neck mobilisation techniques, these versatile tables are suitable for a wide range of physiotherapy and rehabilitation procedures. The "Plus" model features adjustable height armrests which promote a relaxed shoulder girdle and a more comfortable patient during prone treatments. The couch has an excellent height range of 45cm to 98cm and is fitted with a breathing hole and plug as standard.

Models

ST3559 Hydraulic with gas assisted backrest
ST3569 Electric with gas assisted backrest



Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Gas assisted head section adjustable from -25° to +40°
- Adjustable height armrests for prone treatments
- Gas assisted backrest adjustable from 0° to +80°
- Head section and backrest incorporates anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

-EPFS Perimeter Foot Switch

-1WF Single Foot Switch6006 Upper Limb Support

6003 Side Support Rails

9002P Side Support Rail Cushions

6005 Paper Roll Holder

6007 Removable Head/Knee Support



Available Upholstery Colours:

