





Therapy 3 Section Standard Head Couch



A very strong and rigid 3 section couch with a traditionally long head section to suit a wide range of physiotherapy and general examination procedures. The couch can be adjusted to provide support for the patient in prone, supine and seated positions. The couch has an excellent height range of 45cm to 98cm and is fitted with a breathing hole and plug as standard.



ST3557 ST3567



Hydraulic Foot Operation



Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Gas assisted backrest adjustable from 0° to 80°
- Choice of hydraulic or electric height elevation

- Breathing hole and plug included as standard
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- 63cm Wide Easy Clean Upholstery as standard
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

-EPFS Perimeter Foot Switch

Single Foot Switch

Upper Limb Support

Side Support Rails

9002P Side Support Rail Cushions

Paper Roll Holder

Removable Head/Knee Support





Negative head/foot section to -25°

Available Upholstery Colours:

